

Lloyd Davies and Lithotomy positions

Q Can the FLOWTRON® system be used on patients having surgery who are in the Lloyd Davies or Lithotomy positions?

A Yes. Use of *FLOWTRON* calf garments in this patient group is indicated. Patients undergoing surgery in this position should have their limbs frequently checked and repositioned as part of routine surgical care.

The answers to the commonly asked questions may inform and guide but they are not a substitute for holistic patient assessment and clinical judgment. If you cannot find an answer for which you are looking, please register your query at www.act-club.info